

Vectorette PCR of Yeast DNA

Carl Friddle

1) Cut 1-3 μg of clean DNA overnight with 8-10U of blunt cutting enzyme in 20 μl

Most problems come from dirty, uncut DNA. Phenol/glass bead/RNase prepared DNA works well

RsaI, AluI and DraI provide good results.

2) Heat inactivate enzyme

3) Add:

- 3 μl 10x NEBuffer used in digest
- 1 μl annealed anchor bubble
- 1 μl (400U) ligase
- 0.5 μl of 5mM ATP (50 μM ATP final)
- 25.5 μl Water

4) Incubate at 16 C for 9-24 hours.

5) Use 5 μl in 100 μl PCR. Perkin Elmer Ampliwax is recommended for hot start.

- 5 μl of ligation
- 2.5 μl of 20 μM specific primer [M13(-47) for mTn3 library]
- 2.5 μl of 20 μM 224 primer
- 8 μl of 2.5 mM dNTPs
- 10 μl of Taq PCR buffer
- 71 μl Water
- 1 μl Taq DNA polymerase (5U)
- Transfer to Perkin Elmer 9600 Thermal Cycler
 - Denature 92C, 2 minutes
 - 35 Cycles [92C, 20sec; 67C, 30sec; 72C, 45-180sec (>1 min/1 kb)]
 - 72C, 90sec

6) Gel purify 80 μl of PCR product in 1-3% SeaKem GTG, extract with Qiaex (Qiagen), elute with 12 μl of ddWater

7) Sequence 7 μl with Sequenase kit from Amersham.

